

MEOX1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MEOX1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KFS2; MOX1;

Catalog #: A7332

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4222

Isotype: IgG

Swiss Prot: P50221

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a subfamily of non-clustered, diverged, antennapedia-like homeobox-containing genes. The encoded protein may play a role in the molecular signaling network regulating somite development. Alternatively spliced transcript variants encoding different isoforms have been described.

To place an order, please [Click HERE](#).