

RAB5C

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAB5C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RABL; L1880; RAB5L; RAB5CL;

Catalog #: A7342

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5878

Isotype: IgG

Swiss Prot: P51148

Purity: Affinity purification

For research use only.

Background:

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996)

To place an order, please [Click HERE](#).