

ZNF195

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZNF195

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HRF1; ZNFP104;

Catalog #: A7347

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7748

Isotype: IgG

Swiss Prot: O14628

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).