

LZTR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 95kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LZTR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BTBD29; LZTR-1; SWNTS2;

Catalog #: A7350

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8216

Isotype: IgG

Swiss Prot: Q8N653

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the BTB-kelch superfamily. Initially described as a putative transcriptional regulator based on weak homology to members of the basic leucine zipper-like family, the encoded protein subsequently has been shown to localize exclusively to the Golgi network where it may help stabilize the Gogli complex. Deletion of this gene may be associated with DiGeorge syndrome.

To place an order, please [Click HERE](#).