

MAP4K3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 101kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MAP4K3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1;

Catalog #: A7351

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8491

Isotype: IgG

Swiss Prot: Q8IVH8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the mitogen-activated protein kinase kinase kinase family. The encoded protein activates key effectors in cell signalling, among them c-Jun. Alternatively spliced transcripts encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).