

BECN1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BECN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATG6; VPS30; beclin1;

Catalog #: A7353

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8678

Isotype: IgG

Swiss Prot: Q14457

Purity: Affinity purification

For research use only.

Background:

Beclin-1 participates in the regulation of autophagy and has an important role in development, tumorigenesis, and neurodegeneration (Zhong et al., 2009)

To place an order, please [Click HERE](#).