

CER1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CER1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DAND4;

Catalog #: A7354

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9350

Isotype: IgG

Swiss Prot: O95813

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytokine member of the cysteine knot superfamily, characterized by nine conserved cysteines and a cysteine knot region. The cerberus-related cytokines, together with Dan and DRM/Gremlin, represent a group of bone morphogenetic protein (BMP) antagonists that can bind directly to BMPs and inhibit their activity.

To place an order, please [Click HERE](#).