

ARHGAP44

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 89kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ARHGAP44

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RICH2; NPC-A-10;

To place an order, please [Click HERE](#).

Catalog #: A7357

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9912

Isotype: IgG

Swiss Prot: Q17R89

Purity: Affinity purification

For research use only.