

CALCOCO2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CALCOCO2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NDP52;

Catalog #: A7358

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10241

Isotype: IgG

Swiss Prot: Q13137

Purity: Affinity purification

For research use only.

Background:

This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).