

PRDM5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 73kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRDM5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BCS2; PFM2;

Catalog #: A7361

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11107

Isotype: IgG

Swiss Prot: Q9NQX1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcription factor of the PR-domain protein family. It contains a PR-domain and multiple zinc finger motifs. Transcription factors of the PR-domain family are known to be involved in cell differentiation and tumorigenesis.

To place an order, please [Click HERE](#).