

TNFRSF21

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TNFRSF21

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DR6; CD358; BM-018;

Catalog #: A7365

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27242

Isotype: IgG

Swiss Prot: O75509

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the tumor necrosis factor receptor superfamily. The encoded protein activates nuclear factor kappa-B and mitogen-activated protein kinase 8 (also called c-Jun N-terminal kinase 1), and induces cell apoptosis. Through its death domain, the encoded receptor interacts with tumor necrosis factor receptor type 1-associated death domain (TRADD) protein, which is known to mediate signal transduction of tumor necrosis factor receptors. Knockout studies in mice suggest that this gene plays a role in T-helper cell activation, and may be involved in inflammation and immune regulation.

To place an order, please [Click HERE](#).