

MOCS3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MOCS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UBA4;

Catalog #: A7367

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27304

Isotype: IgG

Swiss Prot: O95396

Purity: Affinity purification

For research use only.

Background:

Molybdenum cofactor (MoCo) is necessary for the function of all molybdoenzymes. The protein encoded by this gene adenylates and activates molybdopterin synthase, an enzyme required for biosynthesis of MoCo. This gene contains no introns. A pseudogene of this gene is present on chromosome 14.

To place an order, please [Click HERE](#).