

CARD10

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 116kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CARD10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BIMP1; CARMA3;

Catalog #: A7368

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29775

Isotype: IgG

Swiss Prot: Q9BWT7

Purity: Affinity purification

For research use only.

Background:

The caspase recruitment domain (CARD) is a protein module that consists of 6 or 7 antiparallel alpha helices. It participates in apoptosis signaling through highly specific protein-protein homophilic interactions. Like several other CARD proteins, CARD10 belongs to the membrane-associated guanylate kinase (MAGUK) family and activates NF-kappa-B (NFKB; see MIM 164011) through BCL10 (MIM 603517) (Wang et al., 2001

To place an order, please [Click HERE](#).