

CAMK1G

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CAMK1G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VWS1; CLICK3; CLICKIII; dJ272L16.1;

Catalog #: A7379

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57172

Isotype: IgG

Swiss Prot: Q96NX5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.

To place an order, please [Click HERE](#).