

LGR6

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 104kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LGR6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GPCR; VTS20631;

Catalog #: A7383

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 59352

Isotype: IgG

Swiss Prot: Q9HBX8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the leucine-rich repeat-containing subgroup of the G protein-coupled 7-transmembrane protein superfamily. The encoded protein is a glycoprotein hormone receptor with a large N-terminal extracellular domain that contains leucine-rich repeats important for the formation of a horseshoe-shaped interaction motif for ligand binding. Alternative splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).