

NAA50

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 19kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NAA50

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SAN; MAK3; NAT5; hSAN; NAT13; NAT5P; hNAT5; NAT13P;

To place an order, please [Click HERE](#).

Catalog #: A7387

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80218

Isotype: IgG

Swiss Prot: Q9GZZ1

Purity: Affinity purification

For research use only.