

## CSRNP2

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 60kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CSRNP2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

C12orf2; PPP1R72; TAIP-12; C12orf22; FAM130A1;

**Catalog #:** A7388

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 81566

**Isotype:** IgG

**Swiss Prot:** Q9H175

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene belongs to the CSRNP family of nuclear proteins that share conserved regions, including cysteine- and serine- rich regions, a basic domain, a transcriptional activation domain, and bind the sequence 'AGAGTG', thus have the hallmark of transcription factors. Studies in mice suggest that these genes may have redundant functions. Alternatively spliced transcript variants have been found for this gene.

*To place an order, please [Click HERE](#).*