

RASSF5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RASSF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAPL; Maxp1; NORE1; NORE1A; NORE1B; RASSF3;

Catalog #: A7389

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83593

Isotype: IgG

Swiss Prot: Q8WWW0

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the Ras association domain family. It functions as a tumor suppressor, and is inactivated in a variety of cancers. The encoded protein localizes to centrosomes and microtubules, and associates with the GTP-activated forms of Ras, Rap1, and several other Ras-like small GTPases. The protein regulates lymphocyte adhesion and suppresses cell growth in response to activated Rap1 or Ras. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).