

ANGPT1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human ANGPT1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AGP1; AGPT; ANG1;

Catalog #: A7413

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 284

Isotype: IgG

Swiss Prot: Q15389

Purity: Affinity purification

For research use only.

Background:

Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The protein encoded by this gene is a secreted glycoprotein that activates the receptor by inducing its tyrosine phosphorylation. It plays a critical role in mediating reciprocal interactions between the endothelium and surrounding matrix and mesenchyme and inhibits endothelial permeability. The protein also contributes to blood vessel maturation and stability, and may be involved in early development of the heart. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).