

RIPK1

Reactivity: Human

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RIPK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RIP; RIP1;

To place an order, please [Click HERE](#).

Catalog #: A7414

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8737

Isotype: IgG

Swiss Prot: Q8NFZ0

Purity: Affinity purification

For research use only.