

MAF1

Reactivity: Human Mouse

Tested applications: WB IHC FC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100 FC 1:50 - 1:100

Calculated MW: 29kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human MAF1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7418

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84232

Isotype: IgG

Swiss Prot: Q9H063

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is similar to Maf1, a *Saccharomyces cerevisiae* protein highly conserved in eukaryotic cells. Yeast Maf1 is a negative effector of RNA polymerase III (Pol III). It responds to changes in the cellular environment and represses pol III transcription. Biochemical studies identified the initiation factor TFIIB as a target for Maf1-dependent repression.

To place an order, please [Click HERE](#).