

TEAD3

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TEAD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TEF5; TEAD5; TEF-5; DTEF-1; ETFR-1; TEAD-3;

Catalog #: A7454

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7005

Isotype: IgG

Swiss Prot: Q99594

Purity: Affinity purification

For research use only.

Background:

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon.

To place an order, please [Click HERE](#).