

COQ3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human COQ3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DHHBMT; bA9819.1; DHHBMTASE; UG0215E05;

Catalog #: A7469

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51805

Isotype: IgG

Swiss Prot: Q9NZJ6

Purity: Affinity purification

For research use only.

Background:

Ubiquinone, also known as coenzyme Q, or Q, is a critical component of the electron transport pathways of both eukaryotes and prokaryotes (Jonassen and Clarke, 2000)

To place an order, please [Click HERE](#).