

## ULK4

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:** 142kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human ULK4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FAM7C1; REC01035;

*To place an order, please [Click HERE](#).*

**Catalog #:** A7471

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 54986

**Isotype:** IgG

**Swiss Prot:** Q96C45

**Purity:** Affinity purification

For research use only.