

THAP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human THAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DYT6;

Catalog #: A7472

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55145

Isotype: IgG

Swiss Prot: Q9NVV9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a THAP domain, a conserved DNA-binding domain. This protein colocalizes with the apoptosis response protein PAWR/PAR-4 in promyelocytic leukemia (PML) nuclear bodies, and functions as a proapoptotic factor that links PAWR to PML nuclear bodies. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

To place an order, please [Click HERE](#).