

RRAGC

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RRAGC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GTR2; RAGC; TIB929;

Catalog #: A7479

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64121

Isotype: IgG

Swiss Prot: Q9HB90

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the GTR/RAG GTP-binding protein family. The encoded protein is a monomeric guanine nucleotide-binding protein which forms a heterodimer with RRAGA and RRAGB and is primarily localized to the cytoplasm. The encoded protein promotes intracellular localization of the mTOR complex. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).