

SPAG16

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SPAG16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PF20; WDR29;

Catalog #: A7482

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79582

Isotype: IgG

Swiss Prot: Q8N0X2

Purity: Affinity purification

For research use only.

Background:

Cilia and flagella are comprised of a microtubular backbone, the axoneme, which is organized by the basal body and surrounded by plasma membrane. SPAG16 encodes 2 major proteins that associate with the axoneme of sperm tail and the nucleus of postmeiotic germ cells, respectively (Zhang et al., 2007)

To place an order, please [Click HERE](#).