

CNDP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CNDP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CN1; CPGL2; HsT2308;

Catalog #: A7485

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84735

Isotype: IgG

Swiss Prot: Q96KN2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the M20 metalloprotease family. The encoded protein is specifically expressed in the brain, is a homodimeric dipeptidase which was identified as human carnosinase. This gene contains trinucleotide (CTG) repeat length polymorphism in the coding region.

To place an order, please [Click HERE](#).