

HCCS

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HCCS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLS; CCHL; MCOPS7;

Catalog #: A7490

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3052

Isotype: IgG

Swiss Prot: P53701

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an enzyme that covalently links a heme group to the apoprotein of cytochrome c. Defects in this gene are a cause of microphthalmia syndromic type 7 (MCOPS7). Three transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).