

SP110

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 78kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SP110

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

IPR1; VODI; IFI41; IFI75;

Catalog #: A7492

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3431

Isotype: IgG

Swiss Prot: Q9HB58

Purity: Affinity purification

For research use only.

Background:

The nuclear body is a multiprotein complex that may have a role in the regulation of gene transcription. This gene is a member of the SP100/SP140 family of nuclear body proteins and encodes a leukocyte-specific nuclear body component. The protein can function as an activator of gene transcription and may serve as a nuclear hormone receptor coactivator. In addition, it has been suggested that the protein may play a role in ribosome biogenesis and in the induction of myeloid cell differentiation. Alternative splicing has been observed for this gene and three transcript variants, encoding distinct isoforms, have been identified.

To place an order, please [Click HERE](#).