

RAC3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Catalog #: A7498

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5881

Isotype: IgG

Swiss Prot: P60763

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases.

To place an order, please [Click HERE](#).