

ZNF346

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZNF346

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

JAZ; Zfp346;

Catalog #: A7507

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23567

Isotype: IgG

Swiss Prot: Q9UL40

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a nucleolar, zinc finger protein that preferentially binds to double-stranded (ds) RNA or RNA/DNA hybrids, rather than DNA alone. Mutational studies indicate that the zinc finger domains are not only essential for dsRNA binding, but are also required for its nucleolar localization. The encoded protein may be involved in cell growth and survival.

To place an order, please [Click HERE](#).