

TRIM29

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 66kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRIM29

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATDC;

Catalog #: A7508

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23650

Isotype: IgG

Swiss Prot: Q14134

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the TRIM protein family. It has multiple zinc finger motifs and a leucine zipper motif. It has been proposed to form homo- or heterodimers which are involved in nucleic acid binding. Thus, it may act as a transcriptional regulatory factor involved in carcinogenesis and/or differentiation. It may also function in the suppression of radiosensitivity since it is associated with ataxia telangiectasia phenotype.

To place an order, please [Click HERE](#).