

CSNK1G1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CSNK1G1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CK1gamma1;

Catalog #: A7510

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 53944

Isotype: IgG

Swiss Prot: Q9HCP0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the casein kinase I gene family. The encoded protein is a serine/threonine kinase that phosphorylates acidic proteins, and participates in pathways for cell growth.

To place an order, please [Click HERE](#).