

GOPC

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GOPC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

c

Synonym:

CAL; FIG; PIST; GOPC1; dJ94G16.2;

Catalog #: A7513

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57120

Isotype: IgG

Swiss Prot: Q9HD26

Purity: Affinity purification

For research use only.

Background:

This gene encodes a Golgi protein with a PDZ domain. The PDZ domain is globular and proteins which contain them bind other proteins through short motifs near the C-termini. Mice which are deficient in the orthologous protein have globozoospermia and are infertile. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).