

RAB1B

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAB1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7514

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81876

Isotype: IgG

Swiss Prot: Q9H0U4

Purity: Affinity purification

For research use only.

Background:

Members of the RAB protein family, such as RAB1B, are low molecular mass monomeric GTPases localized on the cytoplasmic surfaces of distinct membrane-bound organelles. RAB1B functions in the early secretory pathway and is essential for vesicle transport between the endoplasmic reticulum (ER) and Golgi (Chen et al., 1997)

To place an order, please [Click HERE](#).