

CASP9

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:10 - 1:100 IF 1:10 - 1:100 IP 1:20 - 1:50

FC 1:20 - 1:50

Calculated MW: 46kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CASP9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

MCH6; APAF3; APAF-3; PPP1R56; ICE-LAP6;

Catalog #: A7523

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 842

Isotype: IgG

Swiss Prot: P55211

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein can undergo autoproteolytic processing and activation by the apoptosome, a protein complex of cytochrome c and the apoptotic peptidase activating factor 1; this step is thought to be one of the earliest in the caspase activation cascade. This protein is thought to play a central role in apoptosis and to be a tumor suppressor. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).