

ZFYVE1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 87kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZFYVE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SR3; DFCP1; TAFF1; ZNFN2A1; PPP1R172;

Catalog #: A7527

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 53349

Isotype: IgG

Swiss Prot: Q9HBF4

Purity: Affinity purification

For research use only.

Background:

The FYVE domain mediates the recruitment of proteins involved in membrane trafficking and cell signaling to phosphatidylinositol 3-phosphate-containing membranes. This protein contains two zinc-binding FYVE domains in tandem and is reported to localize to the Golgi apparatus.

Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).