

## DHRS3

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:1000 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:** 34kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human DHRS3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SDR1; RDH17; Rsdr1; DD83.1; SDR16C1; retSDR1;

**Catalog #:** A7530

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9249

**Isotype:** IgG

**Swiss Prot:** O75911

**Purity:** Affinity purification

For research use only.

**Background:**

Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.

*To place an order, please [Click HERE](#).*