

CBR3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CBR3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

hCBR3; SDR21C2; HEL-S-25;

Catalog #: A7545

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 874

Isotype: IgG

Swiss Prot: O75828

Purity: Affinity purification

For research use only.

Background:

Carbonyl reductase 3 catalyzes the reduction of a large number of biologically and pharmacologically active carbonyl compounds to their corresponding alcohols. The enzyme is classified as a monomeric NADPH-dependent oxidoreductase. CBR3 contains three exons spanning 11.2 kilobases and is closely linked to another carbonyl reductase gene - CBR1.

To place an order, please [Click HERE](#).