

CHN1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CHN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NC; CHN; DURS2; ARHGAP2; RHOGAP2;

Catalog #: A7546

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1123

Isotype: IgG

Swiss Prot: P15882

Purity: Affinity purification

For research use only.

Background:

This gene encodes GTPase-activating protein for ras-related p21-rac and a phorbol ester receptor. It is predominantly expressed in neurons, and plays an important role in neuronal signal-transduction mechanisms. Mutations in this gene are associated with Duane's retraction syndrome 2 (DURS2). Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).