

CSRP2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:21kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CSRP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRP2; LMO5; SmLIM;

Catalog #:A7549

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1466

Isotype:IgG

Swiss Prot:Q16527

Purity:Affinity purification

For research use only.

Background:

CSRP2 is a member of the CSRP family of genes, encoding a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. CRP2 contains two copies of the cysteine-rich amino acid sequence motif (LIM) with putative zinc-binding activity, and may be involved in regulating ordered cell growth. Other genes in the family include CSRP1 and CSRP3. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).