

FCGR2B

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FCGR2B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD32; FCG2; CD32B; FCGR2; IGFR2;

Catalog #: A7554

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2213

Isotype: IgG

Swiss Prot: P31994

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a low affinity receptor for the Fc region of immunoglobulin gamma complexes. The encoded protein is involved in the phagocytosis of immune complexes and in the regulation of antibody production by B-cells. Variations in this gene may increase susceptibility to systemic lupus erythematosus (SLE). Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).