

INHBC

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:38kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human INHBC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IHBC;

Catalog #:A7559

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3626

Isotype:IgG

Swiss Prot:P55103

Purity:Affinity purification

For research use only.

Background:

This gene encodes the beta C chain of inhibin, a member of the TGF-beta superfamily. This subunit forms heterodimers with beta A and beta B subunits. Inhibins and activins, also members of the TGF-beta superfamily, are hormones with opposing actions and are involved in hypothalamic, pituitary, and gonadal hormone secretion, as well as growth and differentiation of various cell types.

To place an order, please [Click HERE](#).