

CCL3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 10kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CCL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MIP1A; SCYA3; G0S19-1; LD78ALPHA; MIP-1-alpha;

Catalog #: A7568

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6348

Isotype: IgG

Swiss Prot: P10147

Purity: Affinity purification

For research use only.

Background:

This locus represents a small inducible cytokine. The encoded protein, also known as macrophage inflammatory protein 1 alpha, plays a role in inflammatory responses through binding to the receptors CCR1, CCR4 and CCR5. Polymorphisms at this locus may be associated with both resistance and susceptibility to infection by human immunodeficiency virus type 1.

To place an order, please [Click HERE](#).