

STK19

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STK19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G11; RP1; D6S60; D6S60E; HLA-RP1;

Catalog #: A7574

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8859

Isotype: IgG

Swiss Prot: P49842

Purity: Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine kinase which localizes predominantly to the nucleus. Its specific function is unknown; it is possible that phosphorylation of this protein is involved in transcriptional regulation. This gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6 and expresses two transcript variants.

To place an order, please [Click HERE](#).