

DIMT1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DIMT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIM1; DIMT1L; HUSSY5; HSA9761;

To place an order, please [Click HERE](#).

Catalog #: A7589

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27292

Isotype: IgG

Swiss Prot: Q9UNQ2

Purity: Affinity purification

For research use only.