

TET3

Reactivity: Human Mouse

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 179kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TET3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

hCG_40738;

Catalog #: A7612

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 200424

Isotype: IgG

Swiss Prot: O43151

Purity: Affinity purification

For research use only.

Background:

Members of the ten-eleven translocation (TET) gene family, including TET3, play a role in the DNA methylation process (Langemeijer et al., 2009)

To place an order, please [Click HERE](#).