

KIF6

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:93kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KIF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C6orf102; dJ137F1.4; dJ188D3.1; dJ1043E3.1;

To place an order, please [Click HERE](#).

Catalog #:A7614

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:221458

Isotype:IgG

Swiss Prot:Q6ZMV9

Purity:Affinity purification

For research use only.