

RNF123

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 149kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RNF123

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KPC1; FP1477;

Catalog #: A7618

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 63891

Isotype: IgG

Swiss Prot: Q5XPI4

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a C-terminal RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions, and an N-terminal SPRY domain. This protein displays E3 ubiquitin ligase activity toward the cyclin-dependent kinase inhibitor 1B which is also known as p27 or KIP1.

To place an order, please [Click HERE](#).